



Quality America, Inc.
We Make Compliance Easy!

IV Access Devices Evaluation

Name:	Occupation/Title:	Dept/Unit:
Today's Date:	Product Name/Number of Times Used:	

Please **circle** the most appropriate answer for each question. Not applicable (N/A) may be used if the question does not apply to this particular product.

agree disagree

1. The safety feature can be activated using a one-handed technique.	1	2	3	4	5	N/A
2. The safety feature does not interfere with normal use of this product.	1	2	3	4	5	N/A
3. Use of this product requires you to use the safety feature.	1	2	3	4	5	N/A
4. This product does not require more time to use than a non-safety device.	1	2	3	4	5	N/A
5. The safety feature works well with a wide variety of hand sizes.	1	2	3	4	5	N/A
6. The device allows for rapid visualization of flashback in the catheter or chamber.	1	2	3	4	5	N/A
7. Use of this product does not increase the number of sticks to the patient.	1	2	3	4	5	N/A
8. The product stops the flow of blood after the needle is removed from the catheter (or after the butterfly is inserted) and just prior to line connections or hep-lock capping.	1	2	3	4	5	N/A
9. A clear and unmistakable change (either audible or visible) occurs when the safety feature is activated.	1	2	3	4	5	N/A
10. The safety feature operates reliably.	1	2	3	4	5	N/A
11. The exposed sharp is blunted or covered after use and prior to disposal.	1	2	3	4	5	N/A
12. The product does not need extensive training to be operated correctly.	1	2	3	4	5	N/A

Of the above questions, which three are the most important to **your** safety when using this product?

Are there other questions which you feel should be asked regarding the safety/utility of this product?

Source: Reprinted with permission of Training for Development of Innovative Control Technology Project. June Fisher, M.D. © June 1993, revised August 1998.

