



Quality America, Inc.
We Make Compliance Easy!

Safety Sharps Evaluation

Name:	Occupation/Title:	Dept/Unit:
Today's Date:	Product Name/Number of Times Used:	

Please **circle** the most appropriate answer for each question. Not applicable (N/A) may be used if the question does not apply to this particular product.

agree disagree

	1	2	3	4	5	N/A
1. I can activate the safety feature with one hand.						
2. I can see the tip of the sharp when I need to (even when the safety feature is activated).						
3. It is impossible NOT to use the safety feature.						
4. This product can be used as quickly as I expected.						
5. The product is easy to handle while wearing gloves.						
6. The device offers a good view of any aspirated fluid.						
7. This product will work with all required syringe and needle sizes.						
8. There is a distinct change (audible or visible) when the safety feature is activated.						
9. The safety feature operates reliably.						
10. The exposed sharp is permanently blunted or covered after use.						
11. The device is just as easy to process after use than our current product.						
12. This product is easy to learn and understand.						
13. The design of the product suggests proper use.						
14. It is almost impossible to skip a crucial step in proper use of the device.						
15. This device provides a better alternative to our current product.						

Comments/Concerns: _____

Adapted from "Training for Development of Innovative Control Technology Project," Trauma Foundation, San Francisco General Hospital, San Francisco, CA.